






Overview: model line **FX-LINE**

FX-LINE – recirculating air filter storage cabinets type 90

focused on the storage of flammable liquids and non flammable hazardous materials

FX-PEGASUS-90 model class with wing doors, door open arrest system, automatic door closing (TSA), one-hand operation (AGT) and recirculating air filter system



FX-CLASSIC-90 model class with wing doors, vertical dividing wall and recirculating air filter system



FX-LINE – recirculating air filter storage cabinets type 30

in addition to the FX-LINE type 90 the stored hazardous materials can be seen at a glance due to a large glass cutout


FX-DISPLAY-30 model class with wing doors with glass openings and recirculating air filter system




CX-LINE – recirculating air filter storage cabinets

focused on the storage of non flammable, slightly aggressive hazardous and toxic materials

CX-CLASSIC-G model class with wing doors and recirculating air filter system




CX-CLASSIC-MultiRisk model class with wing doors and recirculating air filter system



SLX-LINE – recirculating air filter storage cabinets

focused on the storage of acids and alkalis

SLX-CLASSIC model class with wing doors and recirculating air filter system





Flexibility³ with the new recirculating air filter storage cabinets from asecos

Multifunctional

Recirculating air filter storage cabinets from the FX-LINE, CX-LINE and SLX-LINE model groups enable the common storage of the most diverse hazardous materials in one cabinet. Vapours that are harmful to health are extracted and filtered at the same time.



Flexibility³



Flexibility in the choice of the place of installation.



Flexibility when storing different hazardous materials (organic und inorganic) in one cabinet with multi-stage wideband filter.



Flexibility in the choice of the appropriate interior equipment for the recirculating air filter storage cabinet. Select your interior equipment to match the hazardous materials and container sizes stored.

Product lines

FlameFlex

Recirculating air filter storage cabinets from the FX-LINE

The solution for the safe and fire-resistant storage of flammable liquids plus non-flammable hazardous materials. Based on our robust and proven Q-cabinets.

ChemFlex

Recirculating air filter storage cabinets from the CX-LINE

The solution for the safe storage of non-flammable, weakly aggressive hazardous and toxic materials. Based on our high-quality CS-CLASSIC cabinets.

ChemFlex

Recirculating air filter storage cabinets from the SLX-LINE

The solution for the safe storage of non-flammable, weakly aggressive hazardous and toxic materials. Based on our highly resistant SL-CLASSIC cabinets.






Safety at the highest level

Depending on the model, asecos recirculating air filter storage cabinets are able to permanently retain more than 99.998 % of the vapours from hydrocarbons and inorganic hazardous materials up to filter saturation.

The recirculating air filter storage cabinets have been tested and certified in extensive test programmes by infraserv Höchst, an accredited testing body in accordance with GefahrstoffV (German Hazardous Substances Act).



- permanent, monitored technical ventilation of the cabinets, hence no exposure of employees to vapours that are harmful to health.
- elimination of the Ex-zones around the cabinet (when storing flammable liquids)
- elimination of cost-intensive exhaust air systems
- high flexibility in the choice of the place of installation

FX-
CX-
SLX-

Recirculating air filter storage cabinets in daily use



Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel (VDE-tested):

- setting, display and monitoring of the current air recirculation rate
- permanent display and monitoring of the filter saturation (colour display – FlameFlex only)
- display of the time remaining until the next filter change
- reminder alarm for manual filter saturation measurement (inorganic hazardous materials)
- multilingual, pictorial menu guidance
- All alarms are given visually and acoustically. The system additionally has a potential-free alarm contact for connection to a control room.
- device for manually checking the filter saturation



Your safety is our goal!

The recirculating air filter storage cabinets from the FlameFlex and ChemFlex model groups represent the highly efficient combination of regulation-compliant storage of hazardous materials with safe technical ventilation. Further information on the legislation regarding the technical ventilation of safety storage cabinets can be found at www.asecos.com.

The high quality of the asecos recirculating air filter storage cabinets is confirmed by the following certificates and documents:

Certification of the recirculating air filter systems by infraserv Höchst, an accredited testing body in accordance with GefahrstoffV (German Hazardous Substances Act).
Test report no. S084-14

- **fire resistance test of the safety storage cabinet in accordance with EN 14470-1**
- **GS testing of the safety storage cabinet in accordance with the requirements of EN 14470-1 and the laboratory furniture standard EN 14727**
- **VDE approval of the monitoring electronics used**
- **ATEX conformity CE II 3/-G ex ic nA IIB T4 Gc of the recirculating air filter storage cabinet**
- **CE declarations of conformity**

Our certified quality management system in accordance with the international standard ISO 9001 as well as voluntary, permanent external monitoring by a recognised and independent material testing body ensures a constantly high quality standard for all asecos products.





Commissioning and filter change



1 Unpack the device, open the housing front and remove the filter




2 Place the recirculating air filter system on top of the cabinet and connect it to the cabinet ventilation



3 Insert the filter unit, fasten it and close the housing front



4 Plug in the low-temperature IEC power cable and switch on the device



Safety in operation

The filters used for the filtering of inorganic hazardous materials must be checked with all models by means of a manual testing device for filter saturation. This is supported by an electronic reminder function (visual and acoustic). The models from the FX-LINE additionally carry out the automatic electronic detection of hydrocarbons. The filter and exhaust air monitoring is displayed on the touch panel. Visual and acoustic warning signals and a voltage free switch contact provide for maximum safety.



Filter change

Multi-stage wideband filters (e.g. FlameFlex: five-step) with high absorbance capacities are used in all models. An optimised combination of the most diverse, partly impregnated kinds of activated carbon achieves the correspondingly required wideband effect in all devices.

Irrespective of the actual absorbance capacity and the actual saturation of the filter, the exchange of the filter units at least once per year is a major component of our safety concept and of the certification.

The particularly robust design of the filter units permits multiple reprocessing. The asecos system thus works in a particularly environmentally friendly and sustainable manner.



FX
CX
SLX



Overview: model classes **FX-PEGASUS-90**, **FX-CLASSIC-90** and **FX-DISPLAY-30**

Safe and approved storage of flammable liquids and further hazardous substances



Safety for employees

- avoidance of explosive atmospheres
- extraction and filtering of vapours that are harmful to health

Flexibility³

Flexibility³

- flexibility in the choice of the place of installation
- flexibility in the storage of different hazardous materials
- flexibility in the choice of the appropriate interior equipment



Multifunctional

- the FX-LINE recirculating air filter storage cabinets enable the common storage, extraction and filtration of the most diverse hazardous materials in one cabinet

Further features of this model class



Electronic monitoring

- filter saturation by solvent vapours (hydrocarbons)
- air recirculation rate (exhaust air quantity)
- service intervals
- scheduling for the manual filter saturation check



Modern bicolour system

- the new recirculating air filter storage cabinets bring colour to your workplace
- body in modern Anthracite grey
- cabinet doors available in Traffic red, Warning yellow or Light grey



Eco-friendly

- long filter lifetimes
- environmentally and customer-friendly reusable filter system

Possible interior equipment



- **metal-free interior equipment** for the storage of aggressive hazardous materials
- as standard with removable **plastic sump**



- **drawers** with high storage volumes, ideal for storing small containers



- **height-adjustable shelves** in different versions



- **bottom collecting sump** – safe collection of leakages
- usable as storage level in conjunction with perforated insert



Recirculating air filter storage cabinets for the storage of different hazardous substances



- **fire resistance of 90 minutes** (type 90), type-tested in accordance with EN 14470-1



- **fire resistance of 30 minutes** (type 30), type-tested in accordance with EN 14470-1



- **GS approval** of the safety storage cabinet in accordance with EN 14470-1 and the laboratory furniture standard EN 14727



- **excellent products with special manufacturer warranty**
- **extended manufacturer warranty of up to 5 years** in combination with an asecos service tariff

- **VDE approval of the monitoring electronics**



- **effectiveness confirmed** by extensive testing in collaboration with the accredited test center for hazardous materials **infraserV Höchst Co.**
- summary of the test report available on request



FX-
CLASSIC-
SLX



Safe and approved storage of flammable liquids and further hazardous substances

Function / construction:

Cabinet body

- **robust construction and longevity:** triple hinged door, safety elements assembled outside the storage compartment for increased protection against corrosion, scratch- and impact-resistant surface, easy to clean
- **effortless use:** one-hand operation for double-door cabinets (AGT)
- **user-friendly and safe:** door open arrest system (standard), convenient and safe removal or storage of containers
- **always closed:** serial automatic door closing (TSA) releases arrested doors after approx. 60 seconds; acoustic and visual signal before closing
- **no unauthorised use:** doors lockable with cylinder locking (integration in an existing locking system possible) and locking state indicator (red/green)

- **easy transport:** integrated transport base with removable, optional base cover
- **easy alignment:** adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction and filter monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- ATEX compliant: $\text{CE II 3/-G Ex ic nA IIB T4 Gc}$

Recirculating air filter storage cabinet

FX-PEGASUS-90 model FX90.229.120.WDAC

Body colour Anthracite grey (RAL 7016), with wing doors in Traffic red (RAL 3020), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel powder-coated)

Order No. 32905-046-33635

(base cover optional)



Available equipment

- shelves with load capacity of 75 kg
- tray shelves with load capacity of 150 kg
- perforated insert useable as storage area
- bottom collecting sump
- drawers – high storage capacity, optimal use of the available storage capacities
- drawers heavy-duty with load capacity of 100 kg
- metal-free tray shelves and collecting sumps
- storage box
- sump inliner (PE or PP)



Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service life-times and filter recycling



Metal-free interior equipment

Ideal for the storage of flammable, aggressive hazardous substances. Tray shelves with chemical resistant surface lamination and removeable plastic sump, bottom collecting plastic sump.



Easy operation

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



Manual test device

Easy monitoring of the filter saturation when storing inorganic substances



Recirculating air filter storage cabinet

FX-PEGASUS-90 model FX90.229.120.WDAC

Body colour Anthracite grey (RAL 7016), with wing doors in Light grey (RAL 7035), interior equipment with 6 x drawer (sheet steel powder-coated)

Order No. 32905-041-33649

(base cover optional)



Recirculating air filter storage cabinet

FX-PEGASUS-90 model FX90.229.120.WDAC

Body colour Anthracite grey (RAL 7016), with wing doors in Warning yellow (RAL 1004), interior equipment with 4 x drawer (sheet steel powder-coated)

Order No. 32905-040-33647

(base cover optional)

3 door colours are available at no extra cost:



Model

FX90.229.120.WDAC

Information / equipment

without interior equipment

Order No.

32905



Body colour

Anthracite grey RAL 7016
Anthracite grey RAL 7016
Anthracite grey RAL 7016

Door colour

Warning yellow RAL 1004
Light grey RAL 7035
Traffic red RAL 3020

Order No.

040
041
046



Interior equipment packages (can only be ordered in combination with the cabinet)

Material

Order No.

		old code	new code
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated RAL 7035	30003	33635
4x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated RAL 7035	30004	33636
3x shelf, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated RAL 7035	30005	33637
4x shelf, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated RAL 7035	30006	33638
3x tray shelf (V=33.0L), 1x perforated insert, 1x bottom collecting sump (V=55.5L)	sheet steel powder-coated RAL 7035	30007	33640
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	stainless steel	30008	33642
4x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	stainless steel	30009	33643
3x shelf, 1x bottom collecting sump (V=33.0L)	stainless steel	30010	33644
4x shelf, 1x bottom collecting sump (V=33.0L)	stainless steel	30011	33645
3x tray shelf (V=33.0L), 1x bottom collecting sump (V=33.0L)	PP grey	HIGHLIGHT 30013	30013
4x tray shelf (V=33.0L), 1x bottom collecting sump (V=33.0L)	PP grey	HIGHLIGHT 30014	30014
4x drawer (V=30.0L)	sheet steel powder-coated RAL 7035	30015	33647
5x drawer (V=30.0L)	sheet steel powder-coated RAL 7035	30016	33648
6x drawer (V=30.0L)	sheet steel powder-coated RAL 7035	30017	33649
4x drawer (V=30.0L)	stainless steel	30018	33651
5x drawer (V=30.0L)	stainless steel	30019	33652
6x drawer (V=30.0L)	stainless steel	30020	33653

Order example - how to order correctly



32905



040


























































33647



Individual configuration of model **FX90.229.120.WDAC**

 Model	Information / equipment	W x D x H (mm)	Order No.
FX90.229.120.WDAC	without interior equipment	1193 x 615 x 2294	32905

 Body colour	Door colour	Order No.
Anthracite grey RAL 7016	Warning yellow RAL 1004	040 
Anthracite grey RAL 7016	Light grey RAL 7035	041 
Anthracite grey RAL 7016	Traffic red RAL 3020	046 


 Interior equipment	Material	W x D x H (mm)	 1)	 2)	 3)	Order No.	
						old code	new code
 shelf	sheet steel powder-coated RAL 7035	1040 x 432 x 30*	75.00			30021	30021 
	stainless steel	1038 x 425 x 29*	75.00			30024	30024 
 tray shelf	sheet steel powder-coated RAL 7035	1025 x 445 x 71*	150.00	33.00	30.00	30028	30028 
	stainless steel	1025 x 445 x 71*	150.00	33.00	30.00	30064	30064 
 tray shelf, metal-free	PP grey	1002 x 472 x 74*	50.00	33.00	30.00	30027	30027 
 perforated insert	sheet steel powder-coated RAL 7035	1024 x 417 x 60	75.00			5536	5536 
	stainless steel	1024 x 417 x 60	75.00			5773	5773 
 perforated insert, height = 135 mm	sheet steel powder-coated RAL 7035	1024 x 417 x 135	135.00			14084	14084 
 bottom collecting sump	sheet steel powder-coated RAL 7035	1024 x 489 x 86*		33.00	30.00	30022	33634 
	stainless steel	1024 x 489 x 86*	150.00	33.00	30.00	30025	33641 
 bottom collecting sump, V = 55,5 ltr.	sheet steel powder-coated RAL 7035	1017 x 482 x 127*		55.50	50.45	30023	33639 
	stainless steel	1018 x 438 x 120*		55.50	50.00	30026	30026 
 bottom collecting sump, metal-free	PP grey	1002 x 472 x 74*		33.00	30.00	26989	26989 
 drawer	sheet steel powder-coated RAL 7035	844 x 488 x 90*	60.00	30.00	2727	30013	33646 
	stainless steel	844 x 488 x 90*	60.00	30.00	2727	30033	33650 
 drawer (heavy-duty)	sheet steel powder-coated RAL 7035	843 x 487 x 162*	100.00	66.00		30032	30032 
	stainless steel	843 x 487 x 162*	100.00	66.00		30034	30034 
 storage box, left hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360				30029	30029 
 storage box, right hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360				30030	30030 
 sump inliner, for shelf	polyethylene black	1000 x 392 x 50*				7223	7223 
 sump inliner, for drawer	PP grey	834 x 450 x 96*				10054	10054 
 anti-slip rubber mat, for drawer	black	840 x 490 x 1				8699	8699 
 pipe lead-through	sheet steel powder-coated RAL 7035	182 x 104 x 48				30035	30035 
 base cover	sheet steel powder-coated textured RAL 7016	1118 x 35 x 85				30036	30036 
 Q-mover, 2 pieces	steel powder-coated textured RAL 7016	600 x 140 x 130	450.00			30037	30037 
 box for chemicals	PP blue	410 x 300 x 260		18.00	16.36	30038	30038 
 document holder, DIN A5	polystyrene					30793	30793 

* usable storage space W x D (mm)

1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)

Activated carbon replacement filter
with multistage filter system
Order No. 32950 





Technical data


External dimensions W x D x H	mm	1193 x 615 x 2294
Internal dimensions W x D x H	mm	1050 x 522 x 1647
Weight without interior equip.	kg	459
Distributed load	kg/m ²	546
Maximum load	kg	600
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

Transport base

Entry height transport base	mm	90
Entry width transport base	mm	1120

FX90.229.120.WDAC

Top view



FX90.229.120.WDAC

Drawers/number

	9	8	7	6	5	4+1	4	3	2
Distance ⑨ (mm)	150								
Distance ⑧ (mm)	135	150							
Distance ⑦ (mm)	135	170	185						
Distance ⑥ (mm)	135	170	200	250					
Distance ⑤ (mm)	135	170	200	230	250	250			
Distance ④ (mm)	135	170	200	230	295	230	345		
Distance ③ (mm)	135	170	200	230	295	230	360	440	
Distance ② (mm)	170	170	200	230	295	230	360	520	890
Distance ① (mm)	200	200	200	230	295	490	390	520	615

FX90.229.120.WDAC


Shelves/number

	4	3
Distance ④ (mm)	320	
Distance ③ (mm)	350	420
Distance ② (mm)	350	450
Distance ① (mm)	290	380

FX90.229.120.WDAC

Quantity

Quantity	Type of container	Container size
90	Aerosol cans	400 ml
40	Cans	1 litre
40	Laboratory glass bottles	1 litre
40	Laboratory PE-bottles	1 litre
15	Laboratory glass bottles	2.5 litres





Safe and approved storage of flammable liquids and further hazardous substances

Function / construction:

Cabinet body

- **robust construction and longevity:** triple hinged door, safety elements assembled outside the storage compartment for increased protection against corrosion, scratch- and impact-resistant surface, easy to clean
- **effortless use:** one-hand operation for double-door cabinets (AGT)
- **user-friendly and safe:** door open arrest system (standard), convenient and safe removal or storage of containers
- **always closed:** serial automatic door closing (TSA) releases arrested doors after approx. 60 seconds; acoustic and visual signal before closing
- **no unauthorised use:** doors lockable with cylinder locking (integration in an existing locking system possible) and locking state indicator (red/green)

- **easy transport:** integrated transport base with removable, optional base cover
- **easy alignment:** adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction and filter monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- ATEX compliant: **CE** II 3/-G Ex ic nA IIB T4 Gc

Recirculating air filter storage cabinet

FX-PEGASUS-90 model FX90.229.090.WDAC

Body colour Anthracite grey (RAL 7016), with wing doors in Traffic red (RAL 3020), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel powder-coated)

Order No. 32910-046-33622

(base cover optional)



Available equipment

- shelves with load capacity of 75 kg
- tray shelves with load capacity of 150 kg
- perforated insert useable as storage area
- bottom collecting sump
- drawers – high storage capacity, optimal use of the available storage capacities
- drawers heavy-duty with load capacity of 100 kg
- metal-free tray shelves and collecting sumps
- storage box
- sump inliner (PE or PP)



Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service life-times and filter recycling



Easy operation

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



Extraction and filtration

Filtration of solvent vapours (hydrocarbons) as well as inorganic and organic hazardous materials



Manual test device

Easy monitoring of the filter saturation when storing inorganic substances



Recirculating air filter storage cabinet

FX-PEGASUS-90 model FX90.229.090.WDAC

Body colour Anthracite grey (RAL 7016), with wing doors in Light grey (RAL 7035), interior equipment with 4 x drawer (sheet steel powder-coated)

Order No. 32910-041-33627

(base cover optional)



Recirculating air filter storage cabinet

FX-PEGASUS-90 model FX90.229.090.WDAC

Body colour Anthracite grey (RAL 7016), with wing doors in Warning yellow (RAL 1004), interior equipment with 6 x drawer (sheet steel powder-coated)

Order No. 32910-040-33629

(base cover optional)

3 door colours are available at no extra cost:



040 041 046

AVAILABLE from the first quarter of 2016

Order example - how to order correctly



32910 - 040 - 33627

Model	Information / equipment	Order No.
FX90.229.090.WDAC	without interior equipment	32910
Body colour	Door colour	Order No.
Anthracite grey RAL 7016	Warning yellow RAL 1004	040
Anthracite grey RAL 7016	Light grey RAL 7035	041
Anthracite grey RAL 7016	Traffic red RAL 3020	046
Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	old code new code
4x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	30043 33622
3x shelf, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	30044 33623
4x shelf, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	30045 33624
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	stainless steel	30046 33625
4x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	stainless steel	30047 30047
3x shelf, 1x bottom collecting sump (V=22.0L)	stainless steel	30048 30048
4x shelf, 1x bottom collecting sump (V=22.0L)	stainless steel	30049 30049
3x tray shelf (V=22.0L), 1x bottom collecting sump (V=22.0L)	stainless steel	30050 30050
4x tray shelf (V=22.0L), 1x bottom collecting sump (V=22.0L)	PP grey	HIGHLIGHT 30051 30051
4x drawer (V=19.0L)	PP grey	HIGHLIGHT 30052 30052
5x drawer (V=19.0L)	sheet steel powder-coated RAL 7035	30053 33627
6x drawer (V=19.0L)	sheet steel powder-coated RAL 7035	30054 33628
4x drawer (V=19.0L)	sheet steel powder-coated RAL 7035	30055 33629
5x drawer (V=19.0L)	stainless steel	30056 33631
6x drawer (V=19.0L)	stainless steel	30057 33632
	stainless steel	30058 33633



Individual configuration of model **FX90.229.090.WDAC**


Model

FX90.229.090.WDAC

Information / equipment

without interior equipment

W x D x H (mm)

893 x 615 x 2294

Order No.

32910


Body colour

Anthracite grey RAL 7016

Anthracite grey RAL 7016

Anthracite grey RAL 7016

Door colour

Traffic red RAL 3020

Warning yellow RAL 1004

Light grey RAL 7035

Order No.

046

040

041


Interior equipment
Material
W x D x H (mm)

Order No.

old code new code

	shelf	sheet steel powder-coated RAL 7035	740 x 432 x 30*	75.00			30059	30059
		stainless steel	740 x 432 x 30*	75.00			30061	30061
	tray shelf, metal-free	PP grey	702 x 472 x 74*	50.00	23.00	20.91	30063	30063
	perforated insert	sheet steel powder-coated RAL 7035	742 x 417 x 70	75.00			11165	11165
		stainless steel	742 x 417 x 70	75.00			11168	11168
	bottom collecting sump	sheet steel powder-coated RAL 7035	724 x 489 x 83*		22.00	20.00	30060	33621
		stainless steel	724 x 489 x 83*		23.00	20.91	30062	30062
	bottom collecting sump, metal-free	PP grey	702 x 472 x 74*		23.00	20.91	26990	26990
	drawer	sheet steel powder-coated RAL 7035	544 x 488 x 90*	60.00	19.00	1727	30065	33626
		stainless steel	544 x 488 x 90*	60.00	19.00	1727	30067	33630
	drawer (heavy-duty)	sheet steel powder-coated RAL 7035	543 x 487 x 162*	100.00	45.00		30066	30066
		stainless steel	543 x 487 x 162*	100.00	45.00		30068	30068
	storage box, left hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360				30029	30029
	storage box, right hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360				30030	30030
	sump inliner, for drawer	polyethylene white	536 x 480 x 90*		22.00	20.00	27948	27948
	sump inliner, for shelf	PP grey	653 x 406 x 88*				9707	9707
	anti-slip rubber mat, for drawer	black	540 x 490 x 1				15653	15653
	pipe lead-through	sheet steel powder-coated RAL 7035	182 x 104 x 48				30035	30035
	base cover	sheet steel powder-coated textured RAL 7016	818 x 35 x 85				30069	30069
	Q-mover, 2 pieces	steel powder-coated textured RAL 7016	600 x 140 x 130	450.00			30037	30037
	box for chemicals	PP blue	410 x 300 x 260		18.00	16.36	30038	30038
	document holder, DIN A5	polystyrene					30793	30793

* usable storage space W x D (mm)


1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)

Activated carbon replacement filter
with multistage filter system
Order No. 32950







Front view



Side view



Selectional view



Technical data


External dimensions W x D x H	mm	893 x 615 x 2294
Internal dimensions W x D x H	mm	750 x 522 x 1647
Weight without interior equip.	kg	370
Distributed load	kg/m ²	667
Maximum load	kg	600
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

Transport base

Entry height transport base	mm	90
Entry width transport base	mm	820

FX90.229.090.WDAC

Top view



FX90.229.090.WDAC

Drawers/number

	9	8	7	6	5	4+1	4	3	2
Distance ⑨ (mm)	150								
Distance ⑧ (mm)	135	150							
Distance ⑦ (mm)	135	170	185						
Distance ⑥ (mm)	135	170	200	250					
Distance ⑤ (mm)	135	170	200	230	250	250			
Distance ④ (mm)	135	170	200	230	295	230	345		
Distance ③ (mm)	135	170	200	230	295	230	360	440	
Distance ② (mm)	170	170	200	230	295	230	360	520	890
Distance ① (mm)	200	200	200	230	295	490	390	520	615

FX90.229.090.WDAC

Shelves/number

	4	3
Distance ④ (mm)	320	
Distance ③ (mm)	350	420
Distance ② (mm)	350	450
Distance ① (mm)	290	380

Standard storage heights of the drawers

	9	8	7	6	5	4+1	4	3	2
Distance ⑨ (mm)	150								
Distance ⑧ (mm)	135	150							
Distance ⑦ (mm)	135	170	185						
Distance ⑥ (mm)	135	170	200	250					
Distance ⑤ (mm)	135	170	200	230	250	250			
Distance ④ (mm)	135	170	200	230	295	230	345		
Distance ③ (mm)	135	170	200	230	295	230	360	440	
Distance ② (mm)	170	170	200	230	295	230	360	520	890
Distance ① (mm)	200	200	200	230	295	490	390	520	615

Standard storage heights of the shelves

	4	3
Distance ④ (mm)	320	
Distance ③ (mm)	350	420
Distance ② (mm)	350	450
Distance ① (mm)	290	380


FX90.229.090.WDAC

Quantity

49	Aerosol cans	400 ml
23	Cans	1 litre
23	Laboratory glass bottles	1 litre
23	Laboratory PE-bottles	1 litre
9	Laboratory glass bottles	2.5 litres

Storage capacities per drawer/type of container

Type of container	Container size
Aerosol cans	400 ml
Cans	1 litre
Laboratory glass bottles	1 litre
Laboratory PE-bottles	1 litre
Laboratory glass bottles	2.5 litres






Recirculating air filter storage cabinet

FX-PEGASUS-90 model FX90.229.060.WDAC

Body colour Anthracite grey (RAL 7016), with wing door in Traffic red (RAL 3020), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel powder-coated)

Order No. 32915-046-33600

(base cover optional)



Safe and approved storage of flammable liquids and further hazardous substances

Function / construction:

Cabinet body

- **robust construction and longevity:** triple hinged door, safety elements assembled outside the storage compartment for increased protection against corrosion, scratch- and impact-resistant surface, easy to clean
- **user-friendly and safe:** door open arrest system (standard), convenient and safe removal or storage of containers
- **always closed:** serial automatic door closing (TSA) releases arrested door after approx. 60 seconds; acoustic and visual signal before closing
- **no unauthorised use:** door lockable with cylinder locking (integration in an existing locking system possible) and locking state indicator (red/green)

- **easy transport:** integrated transport base with removable, optional base cover
- **easy alignment:** adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction and filter monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- ATEX compliant: II 3/-G Ex ic nA IIB T4 Gc



Manual test device

Easy monitoring of the filter saturation when storing inorganic substances



Easy operation

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



Metal-free interior equipment

Ideal for the storage of flammable, aggressive hazardous substances. Tray shelves with chemical resistant surface lamination and removeable plastic sump, bottom collecting plastic sump.



Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service life-times and filter recycling

Available equipment

- shelves with load capacity of 75 kg
- tray shelves with load capacity of 150 kg
- perforated insert useable as storage area
- bottom collecting sump
- drawers – high storage capacity, optimal use of the available storage capacities
- drawers heavy-duty with load capacity of 100 kg
- metal-free tray shelves and collecting sumps
- storage box
- sump inliner (PE or PP)





Recirculating air filter storage cabinet

FX-PEGASUS-90 model FX90.229.060.WDAC

Body colour Anthracite grey (RAL 7016), with wing door in Warning yellow (RAL 1004), interior equipment with 4 x drawer (sheet steel powder-coated)

Order No. 32915-040-33608

(base cover optional)



Recirculating air filter storage cabinet


FX-PEGASUS-90 model FX90.229.060.WDAC

Body colour Anthracite grey (RAL 7016), with wing door in Light grey (RAL 7035), interior equipment with 6 x drawer (sheet steel powder-coated)

Order No. 32915-041-33610

(base cover optional)

3 door colours are available at no extra cost:



040 041 046

AVAILABLE from the first quarter of 2016

Order example - how to order correctly



Model

Model	Information / equipment	Order No.
FX90.229.060.WDAC	without interior equipment, left hinged	32915
FX90.229.060.WDACR	without interior equipment, right hinged	32920



Body colour

Body colour	Door colour	Order No.
Anthracite grey RAL 7016	Warning yellow RAL 1004	040
Anthracite grey RAL 7016	Light grey RAL 7035	041
Anthracite grey RAL 7016	Traffic red RAL 3020	046



Interior equipment packages (can only be ordered in combination with the cabinet)

	Material	Order No.	
		old code	new code
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035		
4x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	30077	33602
3x shelf, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	30078	33603
4x shelf, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	30079	33604
3x tray shelf (V=12.0L), 1x perforated insert, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated RAL 7035	30080	33606
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	stainless steel	30081	30081
4x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	stainless steel	30082	30082
3x shelf, 1x bottom collecting sump (V=22.0L)	stainless steel	30083	30083
4x shelf, 1x bottom collecting sump (V=22.0L)	stainless steel	30084	30084
3x tray shelf (V=12.0L), 1x perforated insert, 1x bottom collecting sump (V=33.0L)	stainless steel	30085	30085
3x tray shelf (V=11.0L), 1x bottom collecting sump (V=11.0L)	PP grey	HIGHLIGHT 30086	30086
4x tray shelf (V=11.0L), 1x bottom collecting sump (V=11.0L)	PP grey	HIGHLIGHT 30087	30087
4x drawer (V=5.0L)	sheet steel powder-coated RAL 7035	30088	33608
5x drawer (V=5.0L)	sheet steel powder-coated RAL 7035	30089	33609
6x drawer (V=5.0L)	sheet steel powder-coated RAL 7035	30090	33610
3x drawer (V=5.0L), 1x drawer (disposal) (V=30.0L)	sheet steel powder-coated RAL 7035	30091	33612
4x drawer (V=5.0L), 1x drawer (disposal) (V=30.0L)	sheet steel powder-coated RAL 7035	30092	33613
4x drawer (V=5.0L)	stainless steel	30093	33616
5x drawer (V=5.0L)	stainless steel	30094	33617
6x drawer (V=5.0L)	stainless steel	30095	33618
3x drawer (V=5.0L), 1x drawer (disposal) (V=30.0L)	stainless steel	30096	33619
4x drawer (V=5.0L), 1x drawer (disposal) (V=30.0L)	stainless steel	30097	33620



Individual configuration of model **FX90.229.060.WDAC**



Model	Information / equipment	W x D x H (mm)	Order No.
FX90.229.060.WDAC	without interior equipment, left hinged	599 x 615 x 2294	32915
FX90.229.060.WDACR	without interior equipment, right hinged	599 x 615 x 2294	32920



Body colour	Door colour	Order No.
Anthracite grey RAL 7016	Traffic red RAL 3020	046
Anthracite grey RAL 7016	Warning yellow RAL 1004	040
Anthracite grey RAL 7016	Light grey RAL 7035	041



Interior equipment	Material	W x D x H (mm)	1)	2)	3)	Order No.	
						old code	new code
	shelf	sheet steel powder-coated RAL 7035	440 x 432 x 30*	75.00		30098	30098
		stainless steel	440 x 432 x 30*	75.00		30102	30102
	tray shelf	sheet steel powder-coated RAL 7035	427 x 442 x 70*	150.00	12.00	10.91	30099 30099
		stainless steel	427 x 442 x 70*	150.00	12.00	10.91	30103 30103
	tray shelf, metal-free	PP grey	424 x 472 x 67*	50.00	13.00	11.82	30106 30106
	perforated insert	sheet steel powder-coated RAL 7035	420 x 439 x 60	75.00			5555 5555
		stainless steel	425 x 439 x 60	75.00			5776 5776
	bottom collecting sump	sheet steel powder-coated RAL 7035	442 x 507 x 127		22.00	20.00	30100 33594
		stainless steel	418 x 433 x 120*		22.00	20.00	30104 30104
	bottom collecting sump, height = 180 mm	sheet steel powder-coated RAL 7035	409 x 474 x 171*		33.00	30.00	30101 33605
	bottom collecting sump, height = 170 mm	stainless steel	439 x 454 x 170*		33.00	30.00	30105 30105
	bottom collecting sump, metal-free	PP grey	422 x 472 x 67*		12.50	11.36	25424 25424
	drawer	sheet steel powder-coated RAL 7035	345 x 500 x 50*	25.00	5.00	4.55	30108 33607
		stainless steel	345 x 500 x 50*	25.00	5.00	4.55	30110 33614
	drawer (disposal)	sheet steel powder-coated RAL 7035	329 x 462 x 200*	50.00	30.00		30109 30109
		stainless steel	329 x 462 x 150*	50.00	30.00		30111 30111
	storage box, left hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360				30029 30029
	storage box, right hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360				30030 30030
	sump inliner, for drawer disposal	PP grey	318 x 449 x 66*				18472 18472
	sump inliner	PP grey	353 x 406 x 88*				9706 9706
	sump inliner, for drawer disposal	PP black	313 x 413 x 171*				18474 18474
	anti-slip rubber mat, for drawer disposal	black	495 x 340 x 1				8719 8719
	pipe lead-through	sheet steel powder-coated RAL 7035	182 x 104 x 48				30035 30035
	base cover	sheet steel powder-coated textured RAL 7016	524 x 35 x 85				30112 30112
	Q-mover, 2 pieces	steel powder-coated textured RAL 7016	600 x 140 x 130	450.00			30037 30037
	box for chemicals	PP blue	410 x 300 x 260		18.00	16.36	30038 30038
	document holder, DIN A5	polystyrene					30793 30793

* usable storage space W x D (mm)


1) load capacity (kg), 2) retention capacity (litr), 3) maximum container volume in the cabinet (litr)

Activated carbon replacement filter
with multistage filter system
Order No. 32951







Front view



Side view



Selectional view



Technical data


External dimensions W x D x H	mm	599 x 615 x 2294
Internal dimensions W x D x H	mm	450 x 522 x 1647
Weight without interior equip.	kg	287
Distributed load	kg/m ²	915
Maximum load	kg	600
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

Transport base

Entry height transport base	mm	90
Entry width transport base	mm	526

FX90.229.060.WDAC

Top view



FX90.229.060.WDAC TIG#002#

	9	8	7	6	5	4+1	4	3+1	3	2
Distance ⑨ (mm)	145									
Distance ⑧ (mm)	135	145								
Distance ⑦ (mm)	135	170	180							
Distance ⑥ (mm)	135	170	200	240						
Distance ⑤ (mm)	135	170	200	230	240	145				
Distance ④ (mm)	135	170	200	230	295	230	340	275		
Distance ③ (mm)	135	170	200	230	295	230	360	295	435	
Distance ② (mm)	170	170	200	230	295	230	360	295	520	880
Distance ① (mm)	200	200	200	230	295	475	390	475	520	615

FX90.229.060.WDAC Shelves/number


	4	3
Distance ④ (mm)	320	
Distance ③ (mm)	350	420
Distance ② (mm)	350	450
Distance ① (mm)	290	380

FX90.229.060.WDAC Quantity

Quantity	Type of container	Container size
30	Aerosol cans	400 ml
12	Cans	1 litre
12	Laboratory glass bottles	1 litre
12	Laboratory PE-bottles	1 litre
6	Laboratory glass bottles	2.5 litres

Standard storage heights of the shelves

Storage capacities per drawer/type of container





Safe and approved storage of flammable liquids and further hazardous substances

Function / construction:

Cabinet body

- **robust construction and longevity:** triple hinged door, safety elements assembled outside the storage compartment for increased protection against corrosion, scratch- and impact-resistant surface, easy to clean
- **two in one:** 2 compartments (vertical separation) allow the separated storage of different hazardous materials in one cabinet
- **easy handling with comfort:** doors stop in any position, self-closing in the event of fire
- **no unauthorised use:** doors lockable with cylinder locking (integration in an existing locking system possible) and locking state indicator (red/green)

- **easy transport:** integrated transport base with removable, optional base cover
- **easy alignment:** adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction and filter monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- ATEX compliant: $\text{CE II 3/-G Ex ic nA IIB T4 Gc}$

Recirculating air filter storage cabinet FX-CLASSIC-90 model FX90.229.120.MV

Body colour Anthracite grey (RAL 7016), with wing doors in Traffic red (RAL 3020), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel powder-coated), 3 x tray shelf, 1 x bottom collecting sump (PP)

Order No. 33286-046-33656

(base cover optional)

Left compartment



Right compartment




Available equipment

Compartment for the storage of flammable liquids

- shelves with load capacity of 75 kg
- tray shelves with load capacity of 150 kg
- perforated insert useable as storage area
- bottom collecting sump
- drawers – high storage capacity, optimal use of the available storage capacities
- metal-free tray shelves and collecting sumps
- storage box
- sump inliner (PE or PP)

Compartment for the storage of aggressive, non flammable hazardous materials

- metal-free tray shelves and collecting sumps



Easy operation

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



Manual test device

Easy monitoring of the filter saturation when storing inorganic substances



Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service life-times and filter recycling



Separate storage

Dividing wall allows the spatial separated storage of different hazardous materials



Recirculating air filter storage cabinet

FX-CLASSIC-90 model FX90.229.120.MV

Body colour Anthracite grey (RAL 7016), with wing doors in Light grey (RAL 7035), interior equipment with 3 x tray shelf, 1 x bottom collecting sump (PP), 6 x drawer (sheet steel powder-coated)

Order No. 33286-041-33657

(base cover optional)



Recirculating air filter storage cabinet

FX-CLASSIC-90 model FX90.229.120.MV

Body colour Anthracite grey (RAL 7016), with wing doors in Warning yellow (RAL 1004), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (stainless steel), 3 x tray shelf, 1 x bottom collecting sump (PP)

Order No. 33286-040-33544

(base cover optional)

FX-CLASSIC-SLX

3 door colours are available at no extra cost:



040 041 046

AVAILABLE from the first quarter of 2016

Order example - how to order correctly





























































Model	Information / equipment	Order No.
FX90.229.120.MV	without interior equipment, with vertical dividing wall	33286
Body colour	Door colour	Order No.
Anthracite grey RAL 7016	Warning yellow RAL 1004	040
Anthracite grey RAL 7016	Light grey RAL 7035	041
Anthracite grey RAL 7016	Traffic red RAL 3020	046
Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
left: 3x shelf, 1x perforated Insert, 1x bottom collecting sump (V=22.0L) right: 3x tray shelf (V=13.0L), 1x bottom collecting sump (V=11.0L)	sheet steel powder-coated RAL 7035, PP black	old code 33543 new code 33656
left: 3x shelf, 1x perforated Insert, 1x bottom collecting sump (V=22.0L) right: 3x tray shelf (V=13.0L), 1x bottom collecting sump (V=12.5L)	stainless steel, PP grey	33544 33544
left: 6x drawer (V=5.0L) right: 3x tray shelf (V=13.0L), 1x bottom collecting sump (V=11.0L)	PP black, sheet steel powder-coated RAL 7035	33545 33657
left: 6x drawer (V=5.0L) right: 3x tray shelf (V=13.0L), 1x bottom collecting sump (V=11.0L)	PP black, stainless steel	33546 33658



Individual configuration of model **FX90.229.120.MV**

 Model	Information / equipment	W x D x H (mm)	Order No.
FX90.229.120.MV	without interior equipment, with vertical dividing wall	1193 x 615 x 2294	33286

 Body colour	Door colour	Order No.
Anthracite grey RAL 7016	Warning yellow RAL 1004	040 
Anthracite grey RAL 7016	Light grey RAL 7035	041 
Anthracite grey RAL 7016	Traffic red RAL 3020	046 


 Interior equipment	Material	W x D x H (mm)	 1)	 2)	 3)	Order No.	
						old code	new code
	shelf	sheet steel powder-coated RAL 7035	440 x 432 x 30*	75.00		30098	30098 
		stainless steel	440 x 432 x 30*	75.00		30102	30102 
	tray shelf	sheet steel powder-coated RAL 7035	427 x 442 x 70*	150.00	12.00	30099	30099 
		stainless steel	427 x 442 x 70*	150.00	12.00	30103	30103 
	tray shelf, metal-free	PP grey	424 x 472 x 67*	50.00	13.00	30106	30106 
	perforated insert	sheet steel powder-coated RAL 7035	420 x 439 x 60	75.00		5555	5555 
		stainless steel	425 x 439 x 60	75.00		5776	5776 
	bottom collecting sump	sheet steel powder-coated RAL 7035	442 x 507 x 127	22.00	20.00	30100	33594 
		stainless steel	418 x 433 x 120*	22.00	20.00	30104	30104 
	bottom collecting sump, height = 180 mm	sheet steel powder-coated RAL 7035	409 x 474 x 171*	33.00	30.00	30101	33605 
		stainless steel	439 x 454 x 170*	33.00	30.00	30105	30105 
	bottom collecting sump, height = 170 mm	sheet steel powder-coated RAL 7035	422 x 472 x 67*	12.50	11.36	25424	25424 
		stainless steel	422 x 472 x 67*	12.50	11.36	25424	25424 
	drawer	sheet steel powder-coated RAL 7035	345 x 500 x 50*	25.00	5.00	30108	33607 
		stainless steel	345 x 500 x 50*	25.00	5.00	30110	33614 
	drawer (disposal)	sheet steel powder-coated RAL 7035	329 x 462 x 200*	50.00	30.00	30109	30109 
		stainless steel	329 x 462 x 150*	50.00	30.00	30111	30111 
	storage box, left hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360			30029	30029 
	storage box, right hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360			30030	30030 
	sump inliner, for drawer disposal	PP grey	318 x 449 x 66*			18472	18472 
	sump inliner	PP grey	353 x 406 x 88*			9706	9706 
	sump inliner, for drawer disposal	PP black	313 x 413 x 171*			18474	18474 
	anti-slip rubber mat, for drawer disposal	black	495 x 340 x 1			8719	8719 
	pipe lead-through	sheet steel powder-coated RAL 7035	182 x 104 x 48			30035	30035 
	transport base, height = 118 mm	steel powder-coated RAL 7035	1195 x 610 x 118			30171	30171 
	base cover	sheet steel powder-coated textured RAL 7016	1118 x 35 x 85			30036	30036 
	Q-mover, 2 pieces	steel powder-coated textured RAL 7016	600 x 140 x 130	450.00		30037	30037 
	box for chemicals	PP blue	410 x 300 x 260	18.00	16.36	30038	30038 
	document holder, DIN A5	polystyrene				30793	30793 

* usable storage space W x D (mm)

1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)

Activated carbon replacement filter
with multistage filter system
Order No. 32950 




**Technical data**

External dimensions W x D x H	mm	1193 x 615 x 2294
Internal dimensions W x D x H	mm	left: 450 x 522 x 1647 right: 450 x 522 x 1647
Weight without interior equip.	kg	554
Distributed load	kg/m ²	595
Maximum load	kg	600
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

Transport base

Entry height transport base	mm	90
Entry width transport base	mm	1120

FX90.229.120.MV**Top view****FX90.229.120.MV****Drawers/number**

	9	8	7	6	5	4+1	4	3+1	3	2
Distance ⑨ (mm)	145									
Distance ⑧ (mm)	135	145								
Distance ⑦ (mm)	135	170	180							
Distance ⑥ (mm)	135	170	200	240						
Distance ⑤ (mm)	135	170	200	230	240	145				
Distance ④ (mm)	135	170	200	230	295	230	340	275		
Distance ③ (mm)	135	170	200	230	295	230	360	295	435	
Distance ② (mm)	170	170	200	230	265	230	360	295	520	880
Distance ① (mm)	200	200	200	230	295	475	390	475	520	615

FX90.229.120.MV**Shelves/number**

	4	3
Distance ④ (mm)	320	
Distance ③ (mm)	350	420
Distance ② (mm)	350	450
Distance ① (mm)	290	380

FX90.229.120.MV**Quantity**

Quantity	Type of container	Container size
30	Aerosol cans	400 ml
12	Cans	1 litre
12	Laboratory glass bottles	1 litre
12	Laboratory PE-bottles	1 litre
6	Laboratory glass bottles	2.5 litres

Standard storage heights of the drawers (right or left side)


	9	8	7	6	5	4+1	4	3+1	3	2
Distance ⑨ (mm)	145									
Distance ⑧ (mm)	135	145								
Distance ⑦ (mm)	135	170	180							
Distance ⑥ (mm)	135	170	200	240						
Distance ⑤ (mm)	135	170	200	230	240	145				
Distance ④ (mm)	135	170	200	230	295	230	340	275		
Distance ③ (mm)	135	170	200	230	295	230	360	295	435	
Distance ② (mm)	170	170	200	230	265	230	360	295	520	880
Distance ① (mm)	200	200	200	230	295	475	390	475	520	615

Standard storage heights of the shelves (right or left side)

	4	3
Distance ④ (mm)	320	
Distance ③ (mm)	350	420
Distance ② (mm)	350	450
Distance ① (mm)	290	380

Storage capacities per drawer/type of container

Type of container	Container size
Aerosol cans	400 ml
Cans	1 litre
Laboratory glass bottles	1 litre
Laboratory PE-bottles	1 litre
Laboratory glass bottles	2.5 litres






Recirculating air filter storage cabinet

FX-DISPLAY-30 model FX30.229.086.WDFW

Body colour Anthracite grey, with wing door incl. glass cutout in Warning yellow (RAL 1004), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel powder-coated)

Order No. 32946-060-33622 

(base cover optional)



Safe and approved storage of flammable liquids and further hazardous substances


Function / construction:

Cabinet body

- **robust construction and longevity:** triple hinged door, cabinet body made of melamine resin and flame retardant special plates, safety elements assembled outside the storage compartment for increased protection against corrosion, doors with metal front as well as scratch- and impact-resistant surface, easy to clean
- **all containers at a glance:** door with large glass cutout allows a quick and easy overview of the stored containers – even with closed doors
- **easy handling with comfort:** door stops in any position, self-closing in the event of fire
- **no unauthorised use:** door lockable with profile cylinder (integration in an existing locking system possible)

- **easy transport:** integrated transport base with removable, optional base cover or adjusting aids with base cover
- **easy alignment:** adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction and filter monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- ATEX compliant:  II 3/-G Ex ic nA IIB T4 Gc

Available equipment

- shelves with load capacity of 75 kg
- tray shelves with load capacity of 150 kg
- perforated insert useable as storage area
- bottom collecting sump
- drawers – high storage capacity, optimal use of the available storage capacities
- drawers heavy-duty with load capacity of 100 kg
- metal-free tray shelves and collecting sumps
- storage box
- sump inliner (PE or PP)



LARGE STORAGE CAPACITY
—
QUICK GLANCE AT THE STORED ITEMS



Easy operation

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service life-times and filter recycling



Recirculating air filter storage cabinet

FX-DISPLAY-30 model FX30.229.086.WDFW

Body colour Anthracite grey, with wing door in Light grey (RAL 7035), interior equipment with 4 x drawer (sheet steel powder-coated)

Order No. 32946-061-33627

(base cover optional)

Recirculating air filter storage cabinet

FX-DISPLAY-30 model FX30.229.086.WDFW

Body colour Anthracite grey, with wing door in Traffic red (RAL 3020), interior equipment with 6 x drawer (sheet steel powder-coated)

Order No. 32946-066-33629

(base cover optional)

3 door colours are available at no extra cost:



AVAILABLE from the first quarter of 2016

Order example - how to order correctly

















































Model	Information / equipment	Order No.
FX30.229.086.WDFW	without interior equipment, with wing door incl. glass cutout, left hinged	32946
Body colour	Door colour	Order No.
Anthracite grey	Warning yellow RAL 1004	060
Anthracite grey	Light grey RAL 7035	061
Anthracite grey	Traffic red RAL 3020	066
Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	old code 30043 new code 33622
4x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	30044 33623
3x shelf, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	30045 33624
4x shelf, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated RAL 7035	30046 33625
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	stainless steel	30047 30047
4x shelf, 1x perforated insert, 1x bottom collecting sump (V=22.0L)	stainless steel	30048 30048
3x shelf, 1x bottom collecting sump (V=22.0L)	stainless steel	30049 30049
4x shelf, 1x bottom collecting sump (V=22.0L)	stainless steel	30050 30050
3x tray shelf (V=22.0L), 1x bottom collecting sump (V=22.0L)	PP grey	HIGHLIGHT 30051 30051
4x tray shelf (V=22.0L), 1x bottom collecting sump (V=22.0L)	PP grey	HIGHLIGHT 30052 30052
4x drawer (V=19.0L)	sheet steel powder-coated RAL 7035	30053 33627
5x drawer (V=19.0L)	sheet steel powder-coated RAL 7035	30054 33628
6x drawer (V=19.0L)	sheet steel powder-coated RAL 7035	30055 33629
4x drawer (V=19.0L)	stainless steel	30056 33631
5x drawer (V=19.0L)	stainless steel	30057 33632
6x drawer (V=19.0L)	stainless steel	30058 33633



Individual configuration of model **FX30.229.086.WDFW**

 Model	Information / equipment	W x D x H (mm)	Order No.
FX30.229.086.WDFW	without interior equipment, with wing door incl. glass cutout, left hinged	860 x 620 x 2290	32946

 Body colour	Door colour	Order No.
Anthracite grey	Warning yellow RAL 1004	060 
Anthracite grey	Light grey RAL 7035	061 
Anthracite grey	Traffic red RAL 3020	066 


 Interior equipment	Material	W x D x H (mm)	 1)	 2)	 3)	Order No.	
						old code	new code
	shelf	sheet steel powder-coated RAL 7035	740 x 432 x 30*	75.00		30059	30059 
		stainless steel	740 x 432 x 30*	75.00		30061	30061 
	tray shelf, metal-free	PP grey	702 x 472 x 74*	50.00	23.00 20.91	30063	30063 
	perforated insert	sheet steel powder-coated RAL 7035	742 x 417 x 70	75.00		11165	11165 
		stainless steel	742 x 417 x 70	75.00		11168	11168 
	bottom collecting sump	sheet steel powder-coated RAL 7035	724 x 489 x 83*	22.00	20.00	30060	33621 
		stainless steel	724 x 489 x 83*	23.00	20.91	30062	30062 
	bottom collecting sump, metal-free	PP grey	702 x 472 x 74*	23.00	20.91	26990	26990 
	drawer	sheet steel powder-coated RAL 7035	544 x 488 x 90*	60.00	19.00 1727	30065	33626 
		stainless steel	544 x 488 x 90*	60.00	19.00 1727	30067	33630 
	drawer (heavy-duty)	sheet steel powder-coated RAL 7035	543 x 487 x 162*	100.00	45.00	30066	30066 
		stainless steel	543 x 487 x 162*	100.00	45.00	30068	30068 
	storage box, left hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360			30029	30029 
	storage box, right hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360			30030	30030 
	sump inliner, for drawer	polyethylene white	536 x 480 x 90*	22.00	20.00	27948	27948 
	sump inliner, for shelf	PP grey	653 x 406 x 88*			9707	9707 
	anti-slip rubber mat, for drawer	black	540 x 490 x 1			15653	15653 
	pipe lead-through	sheet steel powder-coated RAL 7035	182 x 104 x 48			30035	30035 
	base cover, with adjusting device (available for order only in combination with cabinet)	panel melamine resin-coated RAL 7016	788 x NaN			31010	31010 
	base cover	panel melamine resin-coated RAL 7016				31924	31924 
	Q-mover, 2 pieces	steel powder-coated textured RAL 7016	600 x 140 x 130	450.00		30037	30037 
	box for chemicals	PP blue	410 x 300 x 260	18.00	16.36	30038	30038 

* usable storage space W x D (mm)

1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)

Activated carbon replacement filter
with multistage filter system
Order No. 32951 





Technical data


External dimensions W x D x H	mm	860 x 620 x 2290
Internal dimensions W x D x H	mm	750 x 519 x 1626
Weight without interior equip.	kg	245
Distributed load	kg/m ²	608
Maximum load	kg	600
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

Transport base

Entry height transport base	mm	90
Entry width transport base	mm	788

FX30.229.086.WDFW

Top view



FX30.229.086.WDFW

Drawers/number

	9	8	7	6	5	4+1	4	3	2
Distance ⑨ (mm)	130								
Distance ⑧ (mm)	135	130							
Distance ⑦ (mm)	135	170	160						
Distance ⑥ (mm)	135	170	200	225					
Distance ⑤ (mm)	135	170	200	230	225	225			
Distance ④ (mm)	135	170	200	230	295	230	320		
Distance ③ (mm)	135	170	200	230	295	230	360	420	
Distance ② (mm)	170	170	200	230	295	230	360	520	865
Distance ① (mm)	200	200	200	230	295	490	390	520	615

FX30.229.086.WDFW

Shelves/number

	4	3
Distance ④ (mm)	370	
Distance ③ (mm)	285	400
Distance ② (mm)	315	380
Distance ① (mm)	315	410

Standard storage heights of the drawers

9	8	7	6	5	4+1	4	3	2
---	---	---	---	---	-----	---	---	---

Standard storage heights of the shelves

	4	3
Distance ④ (mm)	370	
Distance ③ (mm)	285	400
Distance ② (mm)	315	380
Distance ① (mm)	315	410


FX30.229.086.WDFW

Quantity

Quantity	Type of container	Container size
49	Aerosol cans	400 ml
23	Cans	1 litre
23	Laboratory glass bottles	1 litre
23	Laboratory PE-bottles	1 litre
9	Laboratory glass bottles	2.5 litres

Storage capacities per drawer/type of container

Quantity	Type of container	Container size
49	Aerosol cans	400 ml
23	Cans	1 litre
23	Laboratory glass bottles	1 litre
23	Laboratory PE-bottles	1 litre
9	Laboratory glass bottles	2.5 litres





Overview: model classes **CX-CLASSIC-G** and **CX-CLASSIC-MultiRisk**

Safe and approved storage of waterpollutants, non flammable liquids and toxics in working areas according to TRGS 510



Multifunctional

- the CX-LINE recirculating air filter storage cabinets enable the common storage, extraction and filtration of the most diverse hazardous materials in one cabinet

Flexibility³

Flexibility³

- flexibility in the choice of the place of installation
- flexibility in the storage of different hazardous materials
- flexibility in the choice of the appropriate interior equipment



Extraction and filtration

- filtration of inorganic and/or weakly aggressive hazardous materials
- recirculating air filter tested and certified by accredited testing body (Infraserv Höchst) in accordance with GefahrstoffV (German Hazardous Substances Act)

Further features of this model class



Safety in operation

- The filters used for the filtering of inorganic hazardous materials must be checked with all models by means of a manual testing device for filter saturation. This is supported by an electronic reminder function (visual and acoustic). The check itself takes place using suitable means of testing (e.g. test tube with pump)



Modern bicolour system


- body in modern Anthracite grey
- cabinet doors available in Pure white or Light grey



Eco-friendly

- long filter lifetimes
- environmentally and customer-friendly reusable filter system

Possible interior equipment



- height-adjustable shelves** in different versions



- pull-out shelves** incl. **sump** made of plastic



- bottom collecting sump** – safe collection of leakages
- usable as storage level in conjunction with perforated insert



Recirculating air filter storage cabinets for the storage, extraction and filtration of non-flammable, slightly aggressive chemicals and toxic substances



- **excellent products with special manufacturer warranty**
- **extended manufacturer warranty of up to five years** in connection with an asecos service tariff

- VDE approval **of the monitoring electronics**



- **effectiveness confirmed** by extensive testing in collaboration with the accredited test center for hazardous materials **infraserv Höchst Co.**
- summary of the test report available on request



FX
CX
SLX



Safe and approved storage of waterpollutants, non flammable liquids and toxics in working areas according to TRGS 510

Function / construction:

Cabinet body

• robust construction and longevity:

robust, high-quality, scratch-resistant outer body with powder-coated textured surface

• no unauthorised use: doors

lockable with cylinder locking (integrated in the turning handle)

• easy alignment: adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- manual check of the filter saturation

Recirculating air filter storage cabinet

CX-CLASSIC-G model CX.229.105.WDFW

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Light grey (RAL 7035), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised)

Order No. 32925-041-30508

Technical data

External dimensions W x D x H	mm	1055 x 520 x 2290
Internal dimensions W x D x H	mm	996 x 493 x 1835
Weight without interior equip.	kg	105
Distributed load	kg/m ²	440
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

CX.229.105.WDFW



Available equipment

- shelves – optionally in galvanised or painted version, height adjustable
- perforated insert useable as storage area
- bottom collecting sump – optionally in galvanised or painted version, liquid-tight
- pull-out shelves with removable plastic sump



Easy operation

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service times and filter recycling



Manual test device

Easy monitoring of the filter saturation when storing inorganic substances



Extraction and filtration

Filtration of solvent vapours (hydrocarbons) as well as inorganic and organic hazardous materials





Recirculating air filter storage cabinet CX-CLASSIC-G model CX.229.105.WDFW

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Pure white (RAL 9010), interior equipment with 4 x pull-out shelf (sheet steel/polyethylene powder-coated/-)

Order No. 32925-050-30516 ⁶

Recirculating air filter storage cabinet CX-CLASSIC-G model CX.229.105.WDFW

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Light grey (RAL 7035), interior equipment with 6 x pull-out shelf (sheet steel/polyethylene powder-coated/-)

Order No. 32925-041-30518 ⁶

2 door colours are available at no extra cost:




















041 050

**AVAILABLE from the first
quarter of 2016**

Order example - how to order correctly



	Model	Information / equipment	Order No.
	CX.229.105.WDFW	without interior equipment, with wing doors incl. glass cutouts	32925
	Body colour	Door colour	Order No.
	Anthracite grey RAL 7016	Light grey RAL 7035	041 
	Anthracite grey RAL 7016	Pure white RAL 9010	050 
	Interior equipment packages (can only be ordered in combination with the cabinet)		Material
	3x shelf, 1x perforated insert, 1x bottom collecting sump [V=33.0L]	sheet steel galvanised	30508 
	3x shelf, 1x bottom collecting sump [V=33.0L]	sheet steel galvanised	30509 
	3x shelf, 1x perforated insert, 1x bottom collecting sump [V=33.0L]	sheet steel powder-coated, sheet steel galvanised	30510 
	3x shelf, 1x bottom collecting sump [V=33.0L]	sheet steel powder-coated, sheet steel galvanised	30511 
	3x shelf, 1x bottom collecting sump with grid [V=33.0L]	sheet steel galvanised, sheet steel painted blue	30512 
	3x shelf, 1x bottom collecting sump with grid [V=33.0L]	sheet steel powder-coated, sheet steel painted blue	30513 
	3x shelf, 1x bottom collecting sump (large) with grid [V=66.0L]	sheet steel galvanised, sheet steel painted blue	30514 
	3x shelf, 1x bottom collecting sump (large) with grid [V=66.0L]	sheet steel powder-coated, sheet steel painted blue	30515 
	4x pull-out shelf [V=33.0L]	sheet steel/polyethylene powder-coated/~ grey/black	30516 
	5x pull-out shelf [V=33.0L]	sheet steel/polyethylene powder-coated/~ grey/black	30517 
	6x pull-out shelf [V=33.0L]	sheet steel/polyethylene powder-coated/~ grey/black	30518 
	Interior equipment	Material	W x D x H (mm)
	shelf	sheet steel galvanised	993 x 492 x 30
	shelf	sheet steel powder-coated	993 x 492 x 30
	perforated insert	sheet steel galvanised	976 x 427 x 41
	bottom collecting sump	sheet steel galvanised	981 x 430 x 85*
	bottom collecting sump with grid	sheet steel coated blue	984 x 424 x 106*
	bottom collecting sump (large) with grid	sheet steel coated blue	984 x 382 x 175*
	pull-out shelf, with PP sump	sheet steel/polypropylene grey/black	898 x 418 x 80*
	box for chemicals	PP blue	410 x 300 x 260

* usable storage space W x D (mm)

1) load capacity [kg], 2) retention capacity [ltr], 3) maximum container volume in the cabinet [ltr]



Safe and approved storage of waterpollutants, non flammable liquids and toxics in working areas according to TRGS 510

Function / construction:

Cabinet body

• robust construction and longevity:

robust, high-quality, scratch-resistant outer body with powder-coated textured surface

• no unauthorised use: doors

lockable with cylinder locking (integrated in the turning handle)

• easy alignment: adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- manual check of the filter saturation

Recirculating air filter storage cabinet

CX-CLASSIC-G model CX.229.081.WDFW

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Light grey (RAL 7035), interior equipment with 4 x pull-out shelf (sheet steel/polyethylene powder-coated/-)

Order No. 32929-041-30548

Technical data

External dimensions W x D x H	mm	810 x 520 x 2290
Internal dimensions W x D x H	mm	751 x 493 x 1835
Weight without interior equip.	kg	93
Distributed load	kg/m²	563
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

CX.229.081.WDFW



Easy operation

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



Manual test device

Easy monitoring of the filter saturation when storing inorganic substances

Available equipment

- shelves – optionally in galvanised or painted version, height adjustable
- perforated insert useable as storage area
- bottom collecting sump – optionally in galvanised or painted version, liquid-tight
- pull-out shelves with removable plastic sump



Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service life-times and filter recycling



Extraction and filtration

Filtration of solvent vapours (hydrocarbons) as well as inorganic and organic hazardous materials



Recirculating air filter storage cabinet CX-CLASSIC-G model CX.229.081.WDFW

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Light grey (RAL 7035), interior equipment with 6 x pull-out shelf (sheet steel/polyethylene powder-coated/-)

Order No. 32929-041-30550

Recirculating air filter storage cabinet CX-CLASSIC-G model CX.229.081.WDFW

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Pure white (RAL 9010), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised)

Order No. 32929-050-30542

2 door colours are available at no extra cost:






























041 050

**AVAILABLE from the first
quarter of 2016**

Order example - how to order correctly



	Model	Information / equipment	Order No.				
	CX.229.081.WDFW	without interior equipment, with wing doors incl. glass cutouts	32929				
	Body colour	Door colour	Order No.				
	Anthracite grey RAL 7016	Light grey RAL 7035	041 				
	Anthracite grey RAL 7016	Pure white RAL 9010	050 				
	Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.				
	3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	sheet steel galvanised	30542 				
	3x shelf, 1x bottom collecting sump (V=33.0L)	sheet steel galvanised	30543 				
	3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated RAL 7035, sheet steel galvanised	30544 				
	3x shelf, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated RAL 7035, sheet steel galvanised	30545 				
	3x shelf, 1x bottom collecting sump with grid (V=33.0L)	sheet steel galvanised, sheet steel painted blue	30546 				
	3x shelf, 1x bottom collecting sump with grid (V=33.0L)	sheet steel powder-coated RAL 7035, sheet steel painted blue	30547 				
	4x pull-out shelf (V=22.0L)	sheet steel/polyethylene powder-coated/~ grey/black	30548 				
	5x pull-out shelf (V=22.0L)	sheet steel/polyethylene powder-coated/~ grey/black	30549 				
	6x pull-out shelf (V=22.0L)	sheet steel/polyethylene powder-coated/~ grey/black	30550 				
	Interior equipment	Material	W x D x H (mm)	 1)	 2)	 3)	Order No.
	shelf	sheet steel galvanised	748 x 492 x 30	60.00			30551 
	shelf	sheet steel powder-coated RAL 7035	748 x 492 x 30	60.00			30552 
	perforated insert	sheet steel galvanised	731 x 412 x 41	60.00			30553 
	bottom collecting sump	sheet steel galvanised	736 x 415 x 112*		33.00	30.00	30554 
	bottom collecting sump with grid	sheet steel coated blue	739 x 411 x 129*		33.00	30.00	30555 
	pull-out shelf, with PP sump	sheet steel/polypropylene grey/black	653 x 418 x 80*	30.00	25.00	22.73	30556 
	box for chemicals	PP blue	410 x 300 x 260		18.00	16.36	30038 
	* usable storage space W x D (mm)						
	1) load capacity (kg), 2) retention capacity (litr), 3) maximum container volume in the cabinet (litr)						
	Accessories	Material					Order No.
	activated carbon replacement filter	activated carbon					32953 



Safe, approved and separated storage of waterpollutant, non flammable liquids and toxics in working areas according to TRGS 510

Function / construction:

Cabinet body

- **robust construction and longevity:** robust, high-quality, scratch-resistant outer body with powder-coated textured surface
- **all in one:** different compartments (multiple separation) allow the storage of different hazardous materials in one cabinet
- **no unauthorised use:** doors lockable with cylinder locking (integrated in the turning handle)

- **easy alignment:** adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- manual check of the filter saturation

Recirculating air filter storage cabinet

CX-CLASSIC-MultiRisk model CX.229.105.MM.4WD

Body colour Anthracite grey (RAL 7016), with wing doors in Pure white (RAL 9010), interior equipment with 4 x shelf, 2 x bottom collecting sump, 4 x perforated insert, 2 x bottom collecting sump (sheet steel galvanised)

Order No. 32933-050-30806



Technical data

External dimensions W x D x H	mm	1055 x 520 x 2290
Internal dimensions W x D x H	mm	all compartments: 475 x 495 x 905
Weight without interior equip.	kg	98
Distributed load	kg/m ²	435
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

CX.229.105.MM.4WD

Available equipment

- shelves – optionally in galvanised or painted version, height adjustable
- perforated insert useable as storage area
- bottom collecting sump – optionally in galvanised or painted version, liquid-tight
- pull-out shelves with removable plastic sump




Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service lifetimes and filter recycling



Extraction and filtration

Filtration of solvent vapours (hydrocarbons) as well as inorganic and organic hazardous materials



MultiRisk marking system

Easy and comfortable marking of the storage compartment with the attached marking system



Recirculating air filter storage cabinet

CX-CLASSIC-MultiRisk model CX.229.105.MM.4WD

Body colour Anthracite grey (RAL 7016), with wing doors in Pure white (RAL 9010), interior equipment with 2 x shelf, 1 x bottom collecting sump, 2 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised), 4 x pull-out shelf (sheet steel/polyethylene powder-coated/-)

Order No. 32933-050-30808

Recirculating air filter storage cabinet

CX-CLASSIC-MultiRisk model CX.229.105.MM.4WD

Body colour Anthracite grey (RAL 7016), with wing doors in Pure white (RAL 9010), interior equipment with 8 x pull-out shelf (sheet steel/polyethylene powder-coated/-)

Order No. 32933-050-31947



050

**AVAILABLE from the first
quarter of 2016**

Order example - how to order correctly



32933



050



30806

Model	Information / equipment	Order No.
CX.229.105.MM.4WD	without interior equipment, with horizontal / vertical dividing wall	32933

Body colour	Door colour	Order No.
Anthracite grey RAL 7016	Pure white RAL 9010	050

Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
4x shelf, 2x bottom collecting sump (V=22.0L), 4x perforated insert, 2x bottom collecting sump (V=22.0L)	sheet steel powder-coated, sheet steel galvanised	30805
4x shelf, 2x bottom collecting sump (V=22.0L), 4x perforated insert, 2x bottom collecting sump (V=22.0L)	sheet steel galvanised	30806
top: 2x shelf, 1x bottom collecting sump (V=22.0L), 2x perforated insert, 1x bottom collecting sump (V=22.0L)	sheet steel powder-coated, sheet steel galvanised	30807
bottom: 4x pull-out shelf (V=15.0L)	RAL 7035, sheet steel/polyethylene powder-coated/- RAL 7035	
top: 2x shelf, 1x bottom collecting sump (V=22.0L), 2x perforated insert, 1x bottom collecting sump (V=22.0L)	sheet steel galvanised RAL 7035, sheet steel/polyethylene powder-coated/- RAL 7035	30808
bottom: 4x pull-out shelf (V=15.0L)		
8x pull-out shelf (V=11.0L)	sheet steel/polyethylene powder-coated/- grey/black	31947

Interior equipment	Material	W x D x H (mm)	1)	2)	3)	Order No.
shelf	sheet steel galvanised		50.00			30838
shelf	sheet steel powder-coated		50.00			30837
perforated insert	sheet steel galvanised	464 x 408 x 41	50.00			30584
bottom collecting sump	sheet steel galvanised			22.00	20.00	30840
bottom collecting sump with grid	sheet steel coated blue			22.00	20.00	30841
pull-out shelf, for cabinets with vertical dividing wall	sheet steel/polypropylene grey/black	378 x 418 x 80*	30.00	15.00	13.64	30839
box for chemicals	PP blue	410 x 300 x 260		18.00	16.36	30038

* usable storage space W x D (mm)

1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)

Accessories	Material	Order No.
activated carbon replacement filter	activated carbon	32952



Overview: model class **SLX-CLASSIC**

Safe and approved storage of aggressive, non flammable hazardous materials in workrooms



Multifunctional

- The SLX-CLASSIC recirculating air filter storage cabinets enable the common storage, extraction and filtration of the most diverse aggressive acids, alkalis and further chemicals in one cabinet

Flexibility³

Flexibility³

- flexibility in the choice of the place of installation
- flexibility in the storage of different hazardous materials
- flexibility in the choice of the appropriate interior equipment



Extraction and filtration

- filtration of inorganic and/or weakly aggressive hazardous materials
- recirculating air filter tested and certified by accredited testing body (Infraserv Höchst) in accordance with GefahrstoffV (German Hazardous Substances Act)

Further features of this model class



Safety in operation

- The filters used for the filtering of inorganic hazardous materials must be checked with all models by means of a manual testing device for filter saturation. This is supported by an electronic reminder function (visual and acoustic). The check itself takes place using suitable means of testing (e.g. test tube with pump)



High quality and durable

- outer body made of epoxy-resin-coated high-grade steel plate and inside made of highly resistant melamine-resin-coated special plates



Eco-friendly

- long filter lifetimes
- environmentally and customer-friendly reusable filter system

Possible interior equipment



- **pull-out shelves** incl. **sump** made of plastic



Recirculating air filter storage cabinets for the storage extraction and filtration of acids and alkalis



- **GS approval** of the cabinet for acids and alkalis in accordance with the laboratory furniture standard EN 14727



- **excellent products with special manufacturer warranty**
- **extended manufacturer warranty of up to five years** in connection with an asecos service tariff

- VDE approval **of the monitoring electronics**



- **effectiveness confirmed** by extensive testing in collaboration with the accredited test center for hazardous materials **infraserv Höchst Co.**
- summary of the test report available on request



FX-
CX-
SLX-



Recirculating air filter storage cabinet SLX-CLASSIC model SLX.230.120.MV

in Light grey (RAL 7035), interior equipment with 12 x pull-out shelf (panel/PP)

Order No. 32936-001-30608



Safe and approved storage of aggressive, non flammable hazardous materials in workrooms

Function / construction:

Cabinet body

- **robust construction and longevity:** cabinet body made of powder-coated sheet steel, safety elements assembled outside the storage compartment for increased protection against corrosion
- **two in one:** two hermetically separated storage compartments allow the isolated storage of acids and bases in one cabinet
- **no unauthorised use:** doors lockable with profile cylinder (integration in an existing locking system possible)

- **easy alignment:** adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- manual check of the filter saturation



Available equipment

- shelf with plastic sump (V=33 ltr.) and grid
- bottom collecting sump PE (V=66 ltr.) with grid
- pull-out shelves with removable plastic sump



Easy operation

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service life-times and filter recycling




Manual test device

Easy monitoring of the filter saturation when storing inorganic substances




Pull-out shelf

Sturdy construction with removable plastic sump

**Technical data****SLX.230.120.MV**

External dimensions W x D x H	mm	1197 x 615 x 2300
Internal dimensions W x D x H	mm	all compartments: 485 x 565 x 1865
Weight without interior equip.	kg	195
Distributed load	kg/m ²	204
Maximum load	kg	200
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

Top viewFX
CX
SLX

001

**AVAILABLE from the first
quarter of 2016**

Order example - how to order correctly



32936



001



30614

**Model**
SLX.230.120.MV**Information / equipment**

without interior equipment, with vertical dividing wall

Order No.

32936

**Body colour**
Light grey RAL 7035**Door colour**

Light grey RAL 7035

Order No.

001

**Interior equipment packages (can only be ordered in combination with the cabinet)**left: 4x pull-out shelf (V=11.0L)
right: 4x pull-out shelf (V=11.0L)
left: 6x pull-out shelf (V=11.0L)
right: 6x pull-out shelf (V=11.0L)**Material**

panel/PP

panel/PP

Order No.

30607

30608

**Interior equipment****Material****W x D x H (mm)**

1)



2)



3)

Order No.pull-out shelf, with PP sump
collecting sump, on shelf with grid
collecting sump, with gridpanel/PP grey
polyethylene
polyethylene407 x 468 x 80*
459 x 555 x 167*
454 x 538 x 314*25.00
70.00
70.0018.00
30.00
66.0016.36
27.27
60.0030609
30610
30611

* usable storage space W x D (mm)

1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)

**Accessories**

activated carbon replacement filter

Material

activated carbon

Order No.

32953



Recirculating air filter storage cabinet
SLX-CLASSIC model SLX.230.060.MH
 in Light grey (RAL 7035), interior equipment with 6 x pull-out shelf
 (panel/PP)
Order No. 32939-001-30615



Safe and approved storage of aggressive, non flammable hazardous materials in workrooms

Function / construction:

Cabinet body

- **robust construction and longevity:** cabinet body made of powder-coated sheet steel, safety elements assembled outside the storage compartment for increased protection against corrosion
- **two in one:** two hermetically separated storage compartments allow the isolated storage of acids and bases in one cabinet
- **no unauthorised use:** doors lockable with profile cylinder (integration in an existing locking system possible)

- **easy alignment:** adjusting aids to compensate for uneven floor

Recirculating air filter system

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- manual check of the filter saturation



Available equipment

- shelf with plastic sump (V=33 ltr.) and grid
- bottom collecting sump PE (V=66 ltr.) with grid
- pull-out shelves with removable plastic sump



Easy operation

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



Reusable filter system

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service times and filter recycling



Manual test device

Easy monitoring of the filter saturation when storing inorganic substances




Pull-out shelf

Sturdy construction with removable plastic sump







Front view



Side view



Selectional view




Technical data

SLX.230.060.MH

External dimensions W x D x H	mm	600 x 615 x 2300
Internal dimensions W x D x H	mm	all compartments: 485 x 565 x 914
Weight without interior equip.	kg	152
Distributed load	kg/m ²	363,00
Maximum load	kg	200
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50

Top view

FX
CX
SLX

001

**AVAILABLE from the first
quarter of 2016**

Order example - how to order correctly



32939



001



30614



Model

Model	Information / equipment	Order No.
SLX.230.060.MH	without interior equipment, with horizontal dividing wall, left hinged	32939
SLX.230.060.MH.R	without interior equipment, with horizontal dividing wall, right hinged	32942



Body colour

Body colour	Door colour	Order No.
Light grey RAL 7035	Light grey RAL 7035	001



Interior equipment packages (can only be ordered in combination with the cabinet)

Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
top: 2x pull-out shelf [V=11.0L] bottom: 2x pull-out shelf [V=11.0L]	panel/PP	30614
top: 3x pull-out shelf [V=11.0L] bottom: 3x pull-out shelf [V=11.0L]	panel/PP	30615



Interior equipment

Interior equipment	Material	W x D x H (mm)	1)	2)	3)	Order No.
pull-out shelf, with PP sump	panel/PP grey	407 x 468 x 80*	25.00	18.00	16.36	30609
collecting sump, on shelf with grid	polyethylene	459 x 555 x 167*	70.00	30.00	27.27	30610
collecting sump, with grid	polyethylene	454 x 538 x 314*	70.00	66.00	60.00	30611

* usable storage space W x D (mm)

1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)



Accessories

Accessories	Material	Order No.
activated carbon replacement filter	activated carbon	32953